

Test Case Planning

Group members	Emily Torres Claire Maye Aidan Isak Sadaf Rasooli Raheema Bahiru Britney Osei
----------------------	--

1. Use Case 1 - <<Symptoms/Side Effects Report>>

Test Case ID	TC01 -
Test Objective	Create a new report
Preconditions	User Logged in Medication is already registered
Test steps	<ol style="list-style-type: none"> 1. Select "Create New Report" 2. Enter Medication name, dosage, side effect description, and date 3. Click "Save" 4. Verify the new report appears
Input Values	Medication name, dosage, description, date
Expected results	A new report is created and saved successfully

Test Case ID	TC02 -
Test Objective	Edit existing report
Preconditions	User signed in, and the report already exists
Test steps	<ol style="list-style-type: none"> 1. Select "Edit Report" 2. Select which report to edit out of the available reports 3. Change information under the report 4. Click finish
Input Values	Strings in report
Expected results	Report information has been changed to hold new edits.

Test Case ID	TC03 -
Test Objective	Deleting an existing report
Preconditions	A report was successfully added

Test steps	<ol style="list-style-type: none">1. Open report menu2. Choose the report to delete and then confirm.3. Check the list of reports again for verification on deletion.
Input Values	None
Expected results	The selected report is removed from the list.

Test Case ID	TC04 -
Test Objective	Ensure that all fields are properly filled when saving a medication report
Preconditions	Open program
Test steps	<ol style="list-style-type: none">1. Select “create report” from the main menu2. Enter medication information, missing dosage, and report text3. Attempt to save, and then close the error box4. Enter the dosage box, but still leave the report text blank5. Attempt to save, and then close the error box
Input Values	Medication name, medication dosage
Expected results	The program will send an error reminding the user to fill all information fields, and preventing them from continuing otherwise